**Ranger 1000 Series instruction Manual**

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---

1.	INTRODUCTION .....	EN-2
1.1	Product Description.....	EN-2
1.2	General Features.....	EN-2
1.3	Definition of Signal Warnings and Symbols.....	EN-3
2.	INSTALLATION .....	EN-5
2.1	Unpacking .....	EN-5
2.2	Installing Components.....	EN-5
2.3	Selecting the Location .....	EN-5
2.4	Leveling the Scale.....	EN-5
2.5	Connecting Power.....	EN-6
	2.5.1 AC Power.....	EN-6
	2.5.2 Battery Power.....	EN-6
2.6	Initial Calibration .....	EN-6
3.	OVERVIEW OF CONTROLS AND FUNCTIONS .....	EN-7
3.1	Display Symbols .....	EN-7
3.2	Controls and Functions.....	EN-7
4.	OPERATIONS.....	EN-8
4.1	Switching the Unit On and Off.....	EN-8
4.2	Manual Tare.....	EN-8
4.3	Zero Operation.....	EN-8
4.4	Changing Weighing Units .....	EN-8
4.5	Accumulation.....	EN-9
5.	SCALE SETTINGS .....	EN-10
5.1	Backlight .....	EN-10
5.2	Zero Tracking Range .....	EN-10
5.3	Filtering .....	EN-10
6.	CALIBRATION .....	EN-11
7.	RESTORE TO FACTORY SETTINGS .....	EN-11
8.	TROUBLESHOOTING.....	EN-12
9.	TECHNICAL DATA.....	EN-13
	WARRANTY .....	EN-16

## 1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the Ranger 1000 Series Weighing Scale. Please read the manual completely before installation and operation.

### 1.1 Product Description

The Ranger 1000 Series is an economical compact bench scale that offers accuracy, durability and portability in an easy-to-use design for light industrial general weighing applications.

### 1.2 General Features

- Capacity models: 3kg, 6kg, 15kg, 30kg with maximum displayed resolution of 1:15,000
- Stainless steel weighing pan
- Durable plastic housing
- 1-window, backlit LCD with prominent 25mm high digits
- Fast display  $\leq$  2 seconds of weight and count results
- kg, g, lb and oz weighing units
- Dedicated TARE key for simple and fast operation
- Accumulation function
- Application modes: Weighing and Accumulation
- 80-hour Internal rechargeable battery with power-saving Auto-shut off

## 1.3 Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

### Signal Words

<b>CAUTION</b>	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided.
<b>Attention</b>	For important information about the product
<b>Note</b>	For useful information about the product

### Warning Symbols

General Hazard

Electrical Shock Hazard

Alternating Current

Direct current

## 1.4 Safety Precautions

**CAUTION:** Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

Verify that the AC Adapter's input voltage range and plug type match the local AC power to be used.

Do not position the scale such that it is difficult to disconnect the power cord from the power receptacle.

Make sure that the power cord does not pose a potential obstacle or tripping hazard.

Operate the equipment only under ambient conditions specified in these instructions.

Do not operate the equipment in hazardous or unstable environments.

This scale is for indoor use only.

Use the scale in dry locations only.

Do not drop loads on the pan.

Use only approved accessories and peripherals.

Disconnect the equipment from the power supply when cleaning.

Service should only be performed by authorized personnel.

## 2. INSTALLATION

### 2.1 Unpacking

Unpack and verify that the following components have been included:

- Scale Unit
- Weighing Platform (plastic pan with stainless steel cover)
- Instruction Manual
- AC Adapter

Save the packaging material. This packaging ensures the ideal protection for the storage or transport of the product.

### 2.2 Installing Components


Position the weighing platform pins into the through-holes on the top housing, then set the platform securely into place. Secure the in-use cover with double-sided tape as needed.

### 2.3 Selecting the Location

Operate the scale on a firm, level surface. Avoid locations with rapid temperature changes, excessive air currents, vibrations, electromagnetic fields, heat or direct sunlight.

### 2.4 Leveling the Scale

Adjust the leveling feet until the bubble is centered in the circle of the level indicator (located on the front panel).

**Attention:** Ensure that the scale is level each time its location is changed.

## 2.5 Connecting Power

### 2.5.1 AC Power

Verify that the intended AC power source matches the AC adapter rating. Connect the supplied AC adapter to the power input receptacle underneath the scale. Plug the AC adapter into a properly grounded power outlet.

### 2.5.2 Battery Power

The battery will begin charging with the AC adapter connected accordingly. An LED indicator below and to the right of the display shows the status of battery charging:

- Green – battery is fully charged
- Yellow – battery is partially charged and charging
- Red – battery is nearly discharged

When AC power is not available, the scale will operate on the internal rechargeable battery. The scale will automatically switch to battery operation if there is a power failure or the power cord is removed. Low battery charge is indicated by the low battery annunciator on the display (the scale will operate for approximately 10 hours more before automatically switching off).

Before using the scale for the first time, the internal rechargeable battery should be fully charged for up to 12 hours. A fully charged battery can operate the scale for approximately 80 hours independent of AC power. The scale can be operated during the charging process. The battery is protected against overcharging and the scale can remain connected to the AC power line.

#### Attention:

- The battery must be recharged every 3 months if the scale is not used for a long time.
- Replace the battery if it does not charge or hold a full charge.
- Dispose of the lead acid battery according to local laws and regulations.

#### CAUTION

**BATTERY IS TO BE REPLACED ONLY BY AN AUTHORIZED SERVICE DEALER. RISK OF EXPLOSION CAN OCCUR IF REPLACED WITH THE WRONG TYPE OR CONNECTED IMPROPERLY.**

## 2.6 Initial Calibration


When the scale is operated for the first time, a Span Calibration is recommended to ensure accurate weighing results. Before performing the calibration, be sure to have the appropriate calibration weights.

Refer to Section 6 for Span Calibration procedures.



### 3. OVERVIEW OF CONTROLS AND FUNCTIONS

#### 3.1 Display Symbols

#### 3.2 Controls and Functions


- Zero the display
- Confirm the parameters modification and move to next setting in User Menu setting

- Toggles to the next weighing unit
- Modify the parameter in User Menu setting

- Inputs the weight of the object on the weighing pan as a Tare value
- Return to the last parameter setting without confirming current modification

- Adds the indicated weight or piece counting value into Accumulation memory
- Long press to enter the User Menu setting
- Skips to the next parameter setting or exits to normal operation

## 4. OPERATIONS

### 4.1 Switching the Unit On and Off

The power switch is located underneath the right-hand side of the scale. Push the switch to the "I" position to turn the scale on, and to the "O" position to turn the scale off. Allow 15-30 minutes for the scale to warm up before use.

**Note:** Make sure the weighing pan is empty before turning the scale on.

### 4.2 Manual Tare

A Tare value in memory is indicated by "NET" on the display.

Place the container on the weighing pan (ex. 100g), then press

**TARE**.

The weight is tared and display will show "0".


A rectangular digital display showing the number "100.00" followed by a small "g" unit indicator.A rectangular digital display showing the word "NET" on the left and "0.0" followed by a small "g" unit indicator on the right.

To clear the Tare value, press **TARE** with the pan empty.

### 4.3 Zero Operation

Center of Zero is indicated by "ZERO" on the display.

Press **ZERO**. "-----" will be displayed momentarily before the display is zeroed.

A rectangular digital display showing the number "100.0" followed by a small "g" unit indicator.A rectangular digital display showing the word "ZERO" on the left and "0.0" followed by a small "g" unit indicator on the right.

### 4.4 Changing Weighing Units

Press **UNITS** successively to toggle through the available weighing units.

## 4.5 Accumulation

Place the item to be weighed on the pan (ex. 300g) then press **M+**. The accumulation entry "ACC 01" is displayed momentarily indicating number of accumulations before showing the total weight and then returns to normal weighing mode.

To perform the next accumulation operation, the scale must be returned to zero.

To display the total accumulation times and weight, press **M+** key when the weight on the pan is less than 10d, then the display will show the accumulation times and total weight successively before returning to normal weighing mode after 4 seconds.

To clear the accumulation data, when the display shows the accumulation data press **No.** key to clear the data.

## 5. SCALE SETTINGS

- In normal weighing mode, long press the **Exit** key to enter the User Menu setting.

In Setup mode:

- o Press **No** to modify the parameters.
- o Press **Yes** to accept the displayed setting and proceed to the next parameter
- o Press **Back** to return to last parameter setting without confirming current modification
- o Press **Exit** to proceed to the next parameter without saving any changes

### 5.1 Backlight

Sets the activation mode of the backlight.

The following settings are available: AU (Auto-on with items greater than 9d placed on the pan or any key is pressed; turns off after 5 seconds of inactivity), on, oFF.

#### Note

Press **Yes** to confirm the backlight setting and move to next setting; or press **Back** to return to last setting.

A rectangular digital display showing the characters 'EL' followed by a space and 'AU'.

### 5.2 Zero Tracking Range

Sets the range in which the zero reading is maintained.


The following settings are available: 0.5d, 1d, 2d, 4d (d = scale division)

A rectangular digital display showing the characters 'R2' followed by a space and '2d'.

### 5.3 Filtering

Sets the level in which the Stable indication turns on; the higher the setting, the faster the stabilization time.

The following settings are available: 0, 1, 2, 3 (levels)


A rectangular digital display showing the characters 'F', a space, '1', a space, 'L', a space, and '1'.

## 6. CALIBRATION

For best results, calibrate the scale at regular intervals. Temperature changes, gravity variations, altitude changes and abuse are few reasons why a scale may need recalibration.

When the scale is ideally positioned for operation, enter calibration and proceed as follows:

Long-press **ZERO** and **UNIT** at the same time to initiate calibration. The required calibration weight flashes on the display (ex. 3000g).

At this time, a different calibration weight value can be selected by pressing **UNIT**. Available selections are 1/3, 2/3 and 100% of full capacity. Once the desired value is displayed, place the corresponding weight on the pan. The display flashes until the actual weight is registered and the calibration process ends (the scale re-starts). At this time, remove the weight on the pan before the display count-down sequence ends.

### Notes:

- The calibration can be performed only when the unit is g, kg or lb.
- The calibration weight value can be changed by pressing **UNIT** key when the value is flickering: 1/3 Max, 2/3 Max or 3/3 Max.

## 7. RESTORE TO FACTORY SETTINGS

During the self-testing period after power on, press Tare and UNIT key at the same time to reset the scale. The scale will now be restored to factory settings.

## 8. MAINTENANCE

### 8.1 Cleaning

**CAUTION:** Disconnect the Ranger 1000 Scale from the power supply before cleaning.

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

**Attention:** Do not use solvents, chemicals, alcohol, ammonia or abrasives to clean the housing or control panel. Do not allow liquids to enter the housing.

### 8.2 TROUBLESHOOTING

The following table lists common problems, possible causes and remedies. Should further problems persist, please contact your local servicing dealer.

Symptom	Possible Causes	Remedy
Scale will not turn on	AC power not connected. Battery discharged.	Connect scale to power.
Battery fails to charge fully	Battery defective or past its useful life.	Replace battery by an authorized service dealer.
Weight reading does not stabilize	Unstable environment. Interference under pan.	Ensure a stable environment. Ensure that the weighing pan is unobstructed and free to move.
Scale does not display accurately	Improper calibration.	Calibrate the scale using proper calibration weights.
E1 error code displayed	EPROM data loss.	Replace the scale.
E2 error code displayed	A/D count is not correct.	Have the load cell replaced.
E3 error code displayed	Display capacity exceeded	Ensure display capacity is not exceeded

### 8.3 Service Information

If the troubleshooting section does not resolve or describe your problem, contact your authorized OHAUS service agent. For service assistance or technical support in the United States call toll-free 1-800-526-0659 between 8:00 AM and 5:00 PM EST. An OHAUS product service specialist will be available to provide assistance. Outside the USA, please visit our web site, [www.ohaus.com](http://www.ohaus.com) to locate the OHAUS office nearest you.

## 9. TECHNICAL DATA

The technical data is valid under the following ambient conditions:

- Indoor use only
- Operating temperature: 0°C (32°F) to 40°C (104°F)
- Humidity: Maximum relative humidity 80 % for temperatures up to 30°C decreasing linearly to 50% relative humidity at 40°C
- Altitude: Up to 2000 m
- Mains supply voltage fluctuations: Up to  $\pm 10\%$  of the nominal voltage
- Installation Category: II
- Pollution Degree: 2

Model	R11P3	R11P6	R11P15	R11P30
Capacity x Readability	3 kg x 0.0002 kg 3000 g x 0.2 g 6 lb x 0.0005 lb 96 x 0.01 oz	6 kg x 0.0005 kg 6000 g x 0.5 g 15 lb x 0.001 lb 240 x 0.02 oz	15 kg x 0.001 kg 15000 g x 1 g 30 lb x 0.002 lb 480 x 0.05 oz	30 kg x 0.002 kg 30000 g x 2 g 60 lb x 0.005 lb 960 x 0.1 oz
Maximum Displayed Resolution	1:15000	1:12000	1:15000	1:15000
Construction	Stainless steel pan, plastic housing			
Weighing Units	kg, g, lb, oz			
Application Modes	Weighing, Accumulation			
Display	1 - window backlit LCD display, 25.4 mm/1" High, 6-digit, 7-segment			
Display Indicators	Stability, Center of Zero, Gross, Tare, Battery status, Units			
Hi-Lo Check Indicators	Display with Alert beeper			
Keyboard	4 Function Membrane switch			
Zero Range	4% of Full Scale Capacity			
Tare Range	Full Capacity by Subtraction			
Stabilization Time	$\leq 2$ seconds			
Operating Temperature	0° to 40°C			
Power	100-240 V ~ 0.3A 50/60Hz input 12 V $\Rightarrow$ 1A output AC Adapter Internal rechargeable sealed lead acid battery			
Battery Life	80 hours continuous use with 12 hour recharge time			
Calibration	Automatic external with kg/g mass, factory calibration recovery			
Shipping Protection	Shipping screw to avoid damage to sensitive components			
Safe Overload Capacity	150% of capacity			
Pan Size	294 x 226 mm / 11.6" x 8.9"			
Scale Dimensions W x H x D	325 x 330.5 x 114 mm / 12.8" x 13" x 4.5"			
Shipping Dimensions W x H x D	440 x 360 x 160 mm / 17.3" x 14.2" x 6.3"			
Net Weight	4.2 kg / 9.3lb			
Shipping Weight	5.3kg / 11.7lb			
Other Features	Auto-Zero Tracking, Accumulation / Totalization			

## Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Marking	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> .

### EC Emissions Note

This device complies with EN55011/CISPR 11 Class B Group 1.

### FCC Note

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

### ISO 9001 Registration

In 1994, Ohaus Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritas Quality International (BVQI), confirming that the Ohaus quality management system is compliant with the ISO 9001 standard's requirements. On May 21, 2009, Ohaus Corporation, USA, was re-registered to the ISO 9001:2008 standard.

### Industry Canada Note

This Class A digital apparatus complies with Canadian ICES-001.

**Disposal**

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

The Batteries Directive 2006/66/EC introduces new requirements from September 2008 on removability of batteries from waste equipment in EU Member States. To comply with this Directive, this device has been designed for safe removal of the batteries at end-of-life by a waste treatment facility.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

For disposal instructions in Europe, refer to [www.ohaus.com/weee](http://www.ohaus.com/weee)

Thank you for your contribution to environmental protection.

## LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period, Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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P/N 30077830 C © 2023 Ohaus Corporation, all rights reserved